

# Safe Therapy for Seniors

## During the COVID-19 Pandemic

Promote safe mobility exercise in senior care facilities – at a distance



**BIODEX**

[www.biodex.com](http://www.biodex.com)

**1-800-224-6339**

Int'l 631-924-9000

# Bring Therapy to the Patient

Reduced staff and social distancing due to the COVID-19 pandemic have raised unprecedented barriers to therapy in senior care. With the help of Biodex devices, staff can keep a physical distance, while keeping therapy on track.

## Mobility Assist™

**Allows safe, independent activity time for patients and residents**

It is important to get patients up and walking to overcome circulation problems, disuse atrophy, and stimulate neural pathways. The Biodex Mobility Assist™ is a motorized stand-assist and walking device that brings patients from a seated to standing position with the protection of a safety harness.

- Minimal staff required to operate
- Patients and residents can safely walk independently
- Enables supervision at a distance when required
- Suitable for patients who require moderate-minimal assistance (BMAT Mobility Levels 2 and 3.)



**[www.biodex.com/mobilityassist](http://www.biodex.com/mobilityassist)**



*Harness safely supports patient during stand assist without risk to therapist*

## medBike®

**Can be easily wheeled to patient's room for cardiac and strength exercise**

Designed for people with physical limitations, medBike® enables both upper and lower movement training in one compact device. Aging and older adults, and post-acute patients not suited to traditional cycling equipment can exercise safely using the medBike.

- Accessible from a chair or wheelchair
- Promotes upper and lower body exercise
- Simple to operate for independent in-room use

\* Distributed exclusively by Biodex in the USA & Canada.



*Easy to understand display*

**[www.biodex.com/medbike](http://www.biodex.com/medbike)**

We're proud to support our partners in healthcare.  
To discuss your specific requirements contact (800) 224-6339, [sales@biodex.com](mailto:sales@biodex.com).

# BIODEX

Biodex Medical Systems, Inc.

20 Ramsey Road, Shirley, New York, 11967-4704, Tel: 800-224-6339 (Int'l 631-924-9000), Fax: 631-924-9241, Email: [info@biodex.com](mailto:info@biodex.com), [www.biodex.com](http://www.biodex.com)